



GEISTER®
The better way to operate™

GEISTER

Aortic Valve
Sizers and
Wire Retractors

GEISTER®



In aortic or mitral valve replacement, assessment of the native annular size and selection of the correspondingly sized prosthetic valve is necessary to achieve maximum hemodynamic performance.

With the dual-ended Memory Aortic Valve Sizer the diameter of the orifice can be determined precisely. The shaft of the instrument is made of Nitinol with thermal shape memory effect. Soft pressure allows the surgeon to shape the sizer to fit the patients geometry. The hollow ends allow a perfect view in the annulus. After heat sterilisation the sizer returns to its original straight shape.

Order information:

03-5475 Memory Aortic Valve Sizers, complete Set of all 8 sizes

consisting of:

03-5477 18/20mm 25cm
03-5478 19/21mm 25cm
03-5479 22/24mm 25cm
03-5480 23/25mm 25cm
03-5481 26/28mm 25cm
03-5482 27/29mm 25cm
03-5483 30/32mm 25cm
03-5484 31/33mm 25cm

With our Wire Retractors with their fenestrated blades you can achieve additional retraction during valvuloplastics. It is also useful for checking the leaflets condition and geometry and for palpating.

Order information:

Atrial Wire Retractor

30-0140	14mm, 24cm
30-0141	17mm, 24cm
30-0142	20mm, 24cm



GEISTER®

The better way to operate™

Geister Medizintechnik GmbH
Föhrenstrasse 2
D-78532 Tuttlingen/Germany
Phone +49 (0) 7461 – 966 240
Telefax +49 (0) 7461 – 966 2422
E-mail info@geister.com
Internet www.geister.com

Distributor